


DENTAL  
WHITENING  
SYSTEMS

ITALDENT



IMPROVES  
MOUTH  
WELLNESS  
WITH A  
BRIGHT SMILE



There's nothing  
that light up health,  
happiness and success  
more than a white  
dazzling smile.

So it's not surprising  
that lots of people choose  
to use dental bleach  
to improve their  
teeth appearance.

# RED BEAM Kit

38% HYDROGEN PEROXIDE

REF RK01

**Red Beam Kit** is a bicomponent system for professional dental bleaching treatment, composed of a gel with 38% of hydrogen peroxide that can be also used with dedicated blue-light or laser-light lamps. The peculiarity of the system is the presence of a red-colour indicator that gradually fade, depending on the oxygen quantity liberated during the bleaching treatment.

Each set of the kit is composed of two syringes, one with high-concentration hydrogen peroxide, and one with a gel. after mixing, you get a highly-reactive gel with 38% of hydrogen peroxide and with a viscosity that doesn't allow its slipping from the dental surface. The presence of potassium ions gives the product also a desensitizing action.

In the kit it's also contained the **Beam Gum Bar-**

**rier Light Curing 3 ml**, a urethane dimethacrylate based light-curing resin with a viscosity that is optimised to be easily applied without smudging. Once polymerized, its particular composition ensures an effective protective barrier for gums. Its colour allows a really precise positioning.

## THE KIT CONTAINS:

2 RedBeam 38% Hydrogen Peroxide Component [ A ]

REF. RB02/1

2 RedBeam 38% Hydrogen Peroxide Component [ B ]

REF. RB02/2

1 Beam Gum Barrier Light Curing 3 ml


REF. GB03

2 Brushes

2 Connectors luer lock

4 Dispensig tips

Instructions



# BLUE BEAM Kit

38% HYDROGEN PEROXIDE

REF BK01

**Blue Beam Kit** is a bicomponent system for professional dental bleaching treatment, based on a gel with 38% of hydrogen peroxide that is chemically activated without the use of blue-light or laser-light lamps.

Each set of the kit is composed of two syringes, one with high-concentration hydrogen peroxide, and one with an activator gel. After the blending, you get a highly-reactive gel with 38% of hydrogen peroxide and with a viscosity that doesn't allow its slipping from the dental surface. The presence of potassium ions gives the product also a desensitizing action.

In the kit it's also contained the **Beam Gum Barrier Light Curing 3 ml**, a urethane dimethacrylate based light-curing resin with a viscosity that is optimised to be easily applied without smudging.

Once polymerized, its particular composition ensures an effective protective barrier for gums. Its colour allows a really precise positioning.

## THE KIT CONTAINS:

2 Blu Beam 38% Hydrogen Peroxide Chemically activated Component [ A ]

REF. BB02/1

2 Blu Beam 38% Hydrogen Peroxide Chemically activated Component [ B ]

REF. BB02/2

1 Beam Gum Barrier light curing 3 ml


REF. GB03

4 Dispensing tips

2 Brushes

2 Connectors luer lock

Instructions



# SHINING Kit

18% E 22% CARBAMIDE PEROXIDE

REF SK01/REF SK02

FOR MAINTENANCE  
AND HOMECARE

**Shining Kit 18%** and **Shining Kit 22%** are products that are based on a sophisticated stabilized carbamide peroxide gel, with pharmaceutical purity grade. The gel shows a consistency and an optimal adhesion, that allow in this way a gradual and uniform release of the active substance. The physiological pH value of  $\approx 6.7$  ensures the enamel's protection. The ma-

terial for the preparation of custom tray is soft and thin, in order to minimize the encumbrance sensation caused by the mask, taking up less interocclusal space and reducing any discomfort to patients.

## EACH KIT CONTAINS:


5 Syringes of mono-component

**REF.** SC001/SC002

5 Dispensing tips

2 2 thermoplastic sheets for trays

Instructions



# EXTRA RED BEAM Kit

MONOCOMPONENT DENTAL BLEACHING SYSTEM BASED ON 38% HYDROGEN PEROXIDE USABLE ALSO WITH LASER-LIGHT LAMP


**REF RK02**

**Extra Red Beam kit** is a system for professional monocomponent dental bleaching treatment, based on a gel with 38% of hydrogen peroxide that is optimized in order to be used with laser-light lamps.

The peculiarity of this system is based on its one-phase composition. The used gel doesn't decompose over time under the high concentration of peroxide, in this way the stability is maintained over

time even without storage at low temperatures.

In the kit it's also contained the Beam Gum Barrier Light Curing 3 ml, a urethane dimethacrylate based light-curing resin for gums' isolation. Both the packaging of monocomponent bleach and of the light-curing dam can be used for 4 applications. The presence of potassium ions gives the product also a desensitizing action.



## THE KIT CONTAINS:

1 Extra Red Beam 38% hydrogen peroxide

**REF. RB01**

1 Syringe Beam Gum Barrier light curing

**REF. GB03**

3 Brushes

4 Dispensing tips

Package leaflet



# EXTRA WHITE BEAM Kit

MONOCOMPONENT DENTAL BLEACHING SYSTEM BASED ON 38% HYDROGEN PEROXIDE USABLE ALSO WITH DEDICATED LAMPS

## REF WK01

**Extra White Beam Kit** Extra White Beam Kit is a system for professional dental bleaching treatment, based on a gel with 38% of hydrogen peroxide that can be used also with dedicated lamps (led-lights, etc..)

The peculiarity of this system is based on its one-phase composition. The used gel doesn't decompose over time under the high concentration of peroxide, in this way the stability is maintained over time even without storage at low temperatures.

In the kit it's also contained the Beam Gum Barrier Light Curing 3 ml, a urethane dimethacrylate based light-curing resin for gums' isolation. Both the packaging of monocomponent bleach and of the light-curing dam can be used for 4 applications.

The presence of potassium ions gives the product also a

desensitizing action.

### THE KIT CONTAINS:

1 Extra White Beam 38% hydrogen peroxide

**REF. WB02**


1 Syringe Beam Gum Barrier light curing

**REF. GB03**

3 Brushes

4 Dispensing tips

Package lealeft



# BEAM GUM BARRIER

LIGHT CURING

REF GB03


**Beam Gum Barrier Light Curing** è un composto foto polimerizzabile a base di uretano dimetacrilato a base di luce-curing con una viscosità che è ottimizzata per essere facilmente applicata con estrema precisione.

Dopo la luce-curing, la sua particolare composizione permette di costruire un'ottima barriera protettiva. Il suo colore consente una reale precisione di posizionamento. La quantità contenuta nella siringa può essere utilizzata

per il trattamento di 5-6 pazienti.

## THE KIT CONTAINS:

- 1 Syringe Beam Gum Barrier light curing 3ml
- 4 Dispensing tips
- Package leaflet





# BEAM DESENS GEL


## DESENSITIZING GEL

REF DEB01

**Beam Desens Gel** is a gel that has specially designed for sensitive teeth. Its composition allows to act on the root causes behind dental sensitivity. The product is particularly recommended in cases of dental sensitivity caused by the whitening dental treatment in study. So it can be used both for professional use and be recommended for the treatments to advise to patients.

### THE KIT CONTAINS:

- 1 Beam Desens Gel
- 3 Brushes
- 3 Dispensing tips
- Package leaflet



# BLUE BEAM ENDO

GEL BASED ON HYDROGEN PEROXIDE , OPTIMIZED FOR ENDODONTIC WHITENING OF NON-VITAL TEETH

REF BE001

**Blue Beam Endo** is a gel based on hydrogen peroxide , that is optimized for endodontic whitening of non-vital teeth. Its particular composition is optimized for the Walking bleach's technique.

### THE KIT CONTAINS:

- 1 Blue Beam Endo
- 3 Dispensing tips
- Package leaflet

